

	Product	Name	Code	Application	Packaging Specification	Image	
Ezi-Con Liquid Transfer Accessories	Linear Filter	Ezi-Con Linear Filter (40µm)	CS-LF-40	For sample filtration; compatible with sterile tube welding.	40µm filter, 5pcs/bag		
		Ezi-Con Linear Filter (100µm)	CS-LF-100	For sample filtration; compatible with sterile tube welding.	100µm filter, 5pcs/bag		
		Ezi-Con Linear Filter (70µm)	CS-LF-70	For sample filtration; compatible with sterile tube welding.	70µm filter, 5pcs/bag		
		Ezi-Con Linear Filter (200µm)	CS-LF-200	For sample filtration; compatible with sterile tube welding.	200µm filter, 5pcs/bag		
	Luer/Spike Connector	Ezi-Con Luer-to-Spike Connector	CS-LSC-01	One end: Spike connector Other end: Female Luer lock connector For transferring liquids into containers equipped with a butterfly port or a pierceable rubber septum.	10pcs/bag		
		Ezi-Con Spike-to-Spike Connector	CS-SSC-01	Both ends are spike connectors, which can be used to connect bags or bottles of reagents with latex interfaces.	10pcs/bag		
		Ezi-Con Double-Spike Sampling Connector	CS-DSC-01	Two Spike Connectors in Parallel, with Heat-Sealed Termination.	10pcs/bag		
		Ezi-Con Luer Adapter Tubing	CS-LAT-01	Both ends are luer male/female interfaces, suitable for luer adapter or sterile connection.	10pcs/bag		
	Dual/Triple Extension Tubing	Ezi-Con Dual Extension Tubing	CS-DET-01	Used for parallel sterile connection of 2 samples or reagents.	5pcs/bag		
		Ezi-Con Triple Extension Tubing	CS-TET-01	Used for parallel sterile connection of 3 samples or reagents.	5pcs/bag		
	Sterile Sampler	Ezi-Con Single Sterile Sampler	CS-SSS-01	The sampling bag has a volume of 8mL and is equipped with an air filter at the end, which can be used for sterile sampling in a C-level environment.	5pcs/bag		
		Ezi-Con Triple Sterile Sampler	CS-TSS-01	Each sampling bag has a volume of 8mL and is equipped with an air filter at the end, which can be used for sterile sampling in a C-level environment.	5pcs/bag		
	Sterile Quantitative Sampling set	Ezi-Con Sterile Quantitative Sampling set	CS-SQS-01	Supplied with a 5mL graduated tube, a three-way stopcock, a syringe, and a luer tubing.	5pcs/bag		
		Ezi-Con Sterile Quantitative Sampling set (L)	CS-SQSL-01	Supplied with a 50mL graduated pipette, an air filter, and a luer tubing.	5pcs/bag		

FlowStori Liquid Transfer Bags&Sets	Sampling-Integrated Fluid Bag	FlowStori Liquid Transfer-Sampling Kit (150mL)	CS-LTK-150	150mL storage bag with an integrated sampler and connecting tubing, designed for liquid storage or final formulation, sampling, and sterile tube welding.	5pcs/bag			
		FlowStori Liquid Transfer-Sampling Kit (600mL)	CS-LTK-600	600mL storage bag with an integrated sampler and connecting tubing, designed for liquid storage or final formulation, sampling, and sterile tube welding.	5pcs/bag			
		FlowStori Liquid Transfer-Sampling Kit (3L)	CS-LTK-3000	3L storage bag with an integrated sampler and connecting tubing, designed for liquid storage or final formulation, sampling, and sterile tube welding.	5pcs/bag			
		FlowStori Liquid Transfer-Sampling Kit (5L)	CS-LTK-5000	5L storage bag with an integrated sampler and connecting tubing, designed for liquid storage or final formulation, sampling, and sterile tube welding.	5pcs/bag			
	Liquid Transfer Bag	FlowStori Liquid Transfer Bag (20mL)	CS-LTB-20	20 mL storage bag, compatible with sterile tube welding, with one male Luer lock connector and one female Luer lock connector. 5pcs/bag	5pcs/bag			
		FlowStori Liquid Transfer Bag (150mL)	CS-LTB-150	150 mL storage bag, compatible with sterile tube welding, with one male Luer lock connector and one female Luer lock connector.	5pcs/bag			
		FlowStori Liquid Transfer Bag (600mL)	CS-LTB-600	600 mL storage bag, compatible with sterile tube welding, with one male Luer lock connector and one female Luer lock connector. 5pcs/bag	5pcs/bag			
		FlowStori Liquid Transfer Bag (1L)	CS-LTB-1000	1L storage bag, compatible with sterile tube welding, with one male Luer lock connector and one female Luer lock connector.	5pcs/bag			
		FlowStori Liquid Transfer Bag (3L)	CS-LTB-3000	3L storage bag, compatible with sterile tube welding, with one male Luer lock connector and one female Luer lock connector. 5pcs/bag	5pcs/bag			
		FlowStori Liquid Transfer Bag (5L)	CS-LTB-5000	5L storage bag, compatible with sterile tube welding, with one male Luer lock connector and one female Luer lock connector.	5pcs/bag			
		FlowStori Liquid Transfer Bag (10L)	CS-LTB-10000	10L storage bag, compatible with sterile tube welding, with one male Luer lock connector and one female Luer lock connector.	5pcs/bag			